

## ANTIOXIDANTS IN FOOD




Figure : (1) PG, (2) THBP, (3) TBHQ, (4) NDGA, (5) BHA, (6) Ionox 100, (7) Octyl gallate, (8) BHT, (9) Lauryl gallate.

Many compounds such as antioxidants, food additives, flavonoids and polyphenolics, toxins, contaminants and drugs are electrochemically active and can readily be measured using HPLC-ECD.

As analytes are identified, both by retention time and voltammetric behavior across the array, compounds misidentification and coelution in complex matrices are minimized even following minimal sample preparation.

### Chromatographic conditions

*Column* : LC-18 (4.6x150mm ; 5 $\mu$ m)-Supelco.

*Mobile Phase* :

A (25mM Sodium acetate and 25mM citric acid/methanol (95/5 v/v).

B (25mM Sodium acetate and 25mM citric acid/methanol acetonitrile(20/40/40 v/v/v)

*Flow Rate* : 1.75mL/min.

*Temperature* : 31°C.

*Injection Volume* : 20 $\mu$ L.